

## Agilent Bioanalyzer Request Form

Name	Date Ordered
Institute and group	
Address	
Email	Phone
Fund	Authorising signature

**Service Required**

RNA	Qualitative Range	Samples per chip
<input type="checkbox"/> RNA 6000 Nano	30-500ng/µl (total)	12
<input type="checkbox"/> RNA 6000 Pico	50-5000pg/µl (total)	11
<input type="checkbox"/> Eukaryotic or <input type="checkbox"/> Prokaryotic or <input type="checkbox"/> Insect		
<input type="checkbox"/> Total RNA or <input type="checkbox"/> messenger RNA		

DNA	DNA Sample Type		
	Sizing Range		
<input type="checkbox"/> DNA 1000	25-1000bp	10-50ng/µl	12
<input type="checkbox"/> High Sensitivity DNA	50-7000bp	50-500pg/µl *	11
<input type="checkbox"/> Genomic DNA (Tape Station)		10-100ng/µl	15 (+ladder)

\*up to 10ng/µl for broad size distribution eg. sheared DNA, NGS libraries

Sample submission - 2µl in appropriate range

Special Instructions

---



---



---

Office use only

Sample receipt checked

Processing

Quality Control Check

Data Transfer

Please include the estimated concentration for each sample:

Sample No.	Sample Name	ULN	
		Nanodrop Reading ng/ul	Well Position
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			

If you have more than 12 samples, please fill in the second table:

Sample No.	Sample Name	ULN	
		Nanodrop Reading ng/ul	Well Position
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			

NB. Sample number depends upon Chip Type